

Stairlift Harness


Designed to offer postural support and maintain a good upright sitting position for those with special needs when using a stair lift.

The full harness consists of a waistbelt, shoulder straps and understrap. Alternatively it is available in separate components as to the users requirements.

This is a permanent fitting harness which is fixed to the bodywork of the stair lift with anchorage brackets. For guidance on fitting it is advisable to contact the stair lift manufacture or supplier.

The stair lift harness or lap belt is sometimes used to support patients when toileting. Please note that you will need to provide your own suitable anchorage points for wall mounting.


If a removable / unfixd harness is preferred please consider Models 19 & 27 in the vehicle section.




Anchorage brackets

Adjusters

Shoulder straps loop onto either side of the waistbelt buckle



Understrap loops over waistbelt buckle



Special Foot & Ankle Strap to hold the feet on the footrest (optional extra)

9yrs to Adult

- STAIR1 - Full harness
- STAIR2 - Belt & Shoulder Straps
- STAIR3 - Belt only

Child to 8yrs

- STAIRJ - Full harness

